IN THE CLAIMS:

- 1. (Cancelled)
- 2. (Cancelled)
- 3. (Currently Amended) A method of evaluating the efficiency of a sterilization process on prion proteins comprising:
 - a) subjecting a sufficient amount of at least one prion protein degradation indicator in
 a container to said sterilization process <u>such that degradation of prion proteins</u>
 <u>occurs</u>; and
 - b) determining the level of degradation of said <u>prion protein degradation</u> indicator, wherein said <u>prion protein degradation</u> indicator is transcribed by a gene selected from the group consisting of SUP35, URE2 and HET-s and is in amyloid form and wherein said level of degradation of said <u>prion protein degradation</u> indicator is indicative of the efficiency of said sterilization process.
- 4. (Cancelled)
- 5. (Previously presented) A method of evaluating the efficiency of a sterilization process on prion proteins, comprising:
 - a) subjecting a sufficient amount of at least one prion protein degradation indicator in
 a container to said sterilization process <u>such that degradation of the prion protein</u>
 <u>occurs</u>, and
 - b) determining the level of degradation of said prion protein degradation indicator,

wherein said <u>prion protein degradation</u> indicator is selected from the group consisting of Sup35p, Ure2p, Het-s protein in amyloid form, and combination thereof and wherein said level of degradation of said <u>prion protein degradation</u> indicator is indicative of the efficiency of said sterilization process.

6. (Previously presented) The method according to c1aim 3, wherein said indicator is a purified naturally occurring form, a recombinant form, a mutant, or a fragment thereof, wherein said indicator is insoluble in non-ionic detergents, resistant to proteases' action, and forms amyloid filaments composed of β-sheets.

7. (Cancelled)

- 8. (Previously presented) The method according to claim 3, wherein step b) is performed by determining a weight or a mass, quantifying radicals, colorimetric variations, radiometry, nephelometry, immuno-enzymatic method, Western blotting, dot blotting, radioimmuno assay, circular dichroism, electron microscopy, fluorescent microscopy, Fourier transform infrared spectroscopy (FTIR), Congo red binding, or proteinase digestion.
- 9. (Previously presented) The method according to claim 3, wherein said sterilization process is performed by autoclaving, chemical exposure, dry heating, low temperature plasma gas, ozone-based exposure, or sterilization techniques using alkylating and/or oxidizing sterilizing agents.
- 10. (Previously presented) The method according to claim 3, wherein said chemical exposure is a vapor or a solution selected from the group consisting of detergent, ethylene oxide, protease, sodium hydroxide, and enzyme.

- 11. (Previously presented) The method of claim 3, wherein said amount of indicator of step a) is between 0.1 ng to 100 g.
- 12. (Previously presented) The method of claim 3, wherein said container is of a material selected from the group consisting of paper, glass, borosilicate, metal, polymer, alloy, and composite.
- 13. (Previously presented) The method according to claim 3, wherein said container is porous, permeable, or semi-permeable.
- 14. (Previously presented) The method of claim 6, wherein said indicator is a purified naturally occurring protein in amyloid form in *Saccharomyces cerevisiae* or *Podospora anserin*.
- 15. (Previously presented) The method according to claim 6, wherein the fragment comprises:
 - a. the first 759 nucleotides of SUP35 counted from the A of the initiation codon encoding the peptidic region,
 - b. the region coding for amino acids 2-114 of Sup35p; or
 - c. the first 639 nucleotides of SUP35 counted from the A of the initiation codon.
- 16. (new) The method according to claim 3 wherein the sterilization process comprises ozone treatment.
- 17. (new) The method according to claim 3 wherein the determining the level of degradation of said prion protein degradation indicator is performed by Western Blot analysis.